

# FERNO<sup>TM</sup> Blood Pressure Cuff

Comfortable, latex-free nylon, velcro cuff

Coiled tube extends to 8 feet

Latex free bladder

Textured brass plated release valve



Large, high contrast, easy-to-read 6" dial

Wall mount with 110 degree rotation

## FERNO<sup>TM</sup> Blood Pressure Cuff

If you're looking for an easy-to-use, durable, and reliable tool that accurately reads blood pressures in your facility, Spiracle Technology's wall-mounted blood pressure cuff is the perfect fit. Spiracle Technology's wall-mounted blood pressure cuff is resilient; it withstands the harsh use common in pre-hospital settings. This attractive and affordable blood pressure cuff has a maximum reading of 300 mmHg and is easy-to-read; the 6" dial rotates 110 degrees for easy viewing. In addition to a comfortable, latex-free, nylon Velcro cuff, Spiracle Technology's wall-mounted blood pressure cuff also includes a rear bracket compartment which provides plenty of storage space.

### Description

- Wall mounted
- Large easy to read 6" (150mm) dial
- 300mm Hg No pin Stop Mechanism
- Rear bracket storage compartment area
- 110 degree rotation for easy viewing
- Nylon velcro cuff
- Latex free bladder
- plated brass air release valve
- Coil tubing extends to 8 feet
- Textured deflation valve
- Highly durable
- Compliant to EN 1060
- CE marked

### How to Order

**Model T1046019**  
Aneroid Sphygmomanometer  
(Blood Pressure Cuff)  
Wall mounted type



 **Spiracle Technology**  
A FERNO Group Company

10601 Calle Lee #190  
Los Alamitos, CA 90720-2524 U.S.A.  
Phone: (714) 418-1091 • Fax: (714) 418-1095  
e-mail: info@spiracle.com • www.spiracle.com