

# Immobilization Solutions



Ferno's complete line of **immobilization solutions** provide **fast and safe spinal immobilization** with almost any backboard or scoop stretcher. Ferno immobilization solutions fit both adult and pediatric patients.



**Fastrap Quick Restraint System** quickly secures patients to backboards and other devices



**WizLoc Cervical Collar** is adjustable for patients from age two through adult



**Universal Head Immobilizer** is lightweight, easy to apply and adaptable to any size patient

## Fastrap™ Quick Restraint System



Quickly secure patients to backboards and other devices with the Ferno Fastrap™ Quick Restraint System. This 10-point, six-strap system is sequentially color-coded (black, red, yellow, green and gray) and includes nylon webbing with hook-and-loop closures at the ends of the straps for fast, easy application.

- Size adjustable from pediatric to adult
- Easy storage and rapid application
- Compatible with almost any backboard or scoop stretcher

## Universal Head Immobilizer



The Universal Head Immobilizer is lightweight and easy to apply. The baseplate fits almost any backboard and the head supports can be inverted for application on a Scoop Stretcher.

- Head straps conform to the patient's forehead
- Large ear holes allow easy patient monitoring
- Minimal interference with X-ray, MRI, or CT procedures
- PVC Plastic-coated, closed-cell foam head supports are easy to clean, will not absorb blood or bodily fluids, and prevent bacterial growth
- Meets FMVSS 302 flammability rating

## WizLoc Cervical Collar



The WizLoc Cervical Collar has a universal design that fits a wide range of patient sizes from age two through adult. The adjustable, 3-position occipital support with keyhole lock-pin and self-centering mandibular support provide on-patient sizing and adjustment for increased accuracy of application and greater patient comfort and safety.

- Minimal interference with X-ray or MRI procedures
- Can be decontaminated and is reusable up to 10 times
- Available in high visibility neon pink or OD green

Learn more at [www.fernoems.com/immobilization](http://www.fernoems.com/immobilization)