

# FERNO Compact 1 & 2 Carry Chairs

Ideal for use in confined areas



Compact 1 Carry Chair



Compact 2 Carry Chair

FOR MORE DETAILED  
INFORMATION SEE OVER

**FERNO** – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Carry Chairs** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

## Specifications

	Height (folded)	Width (folded)	Depth (folded)	Weight	Load Limit
Compact 1	935mm (765mm)	435mm (435mm)	430 mm (235mm)	9.1 kg	200kg/31st/440lbs
Compact 2	935mm (775mm)	445mm (445mm)	420 mm (285mm)	9.4 kg	200kg/31st/440lbs

call +44 (0) 1274 851 999  
 email sales@ferno.co.uk  
 www ferno.co.uk



# FERNO Compact 1 & 2 Carry Chairs

The Compact Carry Chair range of chairs are durable, robust and are manufactured from stainless steel. Their one piece back and seat flame retardant, vinyl coated nylon covers allows the patient to adopt a comfortable position during transfer.

The Compact range design incorporates safety rings which slide over the folding mechanism.

- Seat covers can simply be wiped clean for maximum hygiene.
- Folds neatly into a compact size for easy stowage.
- 7" wheels for ease of manoeuvring.
- Maximum Load 200kg.
- Easily adjustable chest and ankle restraints are quickly applied
- Unique style lifting loops help the user to adopt a secure grip
- Large wheels and an axle bar to provide stability over uneven ground

The Compact 1 proven chair is invaluable when transferring patients through corridors and lifts due to it's narrow design and large rear wheels.

The Compact 2 Carry Chair has been designed with ergonomic front grip handles for improved posture when carrying patients down stairs.



#### Part Nos:

071481302/S	Compact 1 Carry Chair (black)
071481401/S	Compact 2 Carry Chair (black)
071481401/BIO	Compact 2 Carry Chair (black) with Biosafe straps

#### Accessories

038299401	Dust cover for Compact 1, & 2 (black)
081931400	Wall Bracket for Compact 1 & 2 chairs
031350600	Carry Case For Compact 1/2 Chair

FERNO Compact 1 & 2 Carry Chairs FPC002 10.14

call +44 (0) 1274 851 999  
email [sales@ferno.co.uk](mailto:sales@ferno.co.uk)  
www [ferno.co.uk](http://ferno.co.uk)



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.