

# Mobyle Chair MK1 & MK4



## **PERSONAL ONSITE DEMO**

These practical, comfortable carrying chairs serve the demanding requirements of the emergency services and are widely used by professionals in other areas, such as Nursing Homes, Hospitals, Clinics, Sports/Leisure Establishments and various types of industry.

Ideally suited for carrying patients up and down stairs the wheels are set inside the frame for minimum width to allow use in confined spaces.

#### Standard features include:

A chest restraint with an Easi-Clip buckle and the raised seat panel enhance patient comfort and security.

Manufactured from steel with powder coated finish. The seat panels and back are made in wipe clean, vinyl coated, flame retardant nylon and may be removed for cleaning. The chairs fold to a slim size for storage.



Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire BD19 4TZ England.









Available in a range of two models designed to offer a choice of features to meet the demands of differing environments.

The Mk1 has two wheels and attendant lifting loops.

The Mk4 has four wheels, a raised foot rest and lifting loops, plus a higher back panel.

### **TECHNICAL DATA - Mobyle Chair**

#### Mobyle MK1 Mobyle MK4 Height 91cm Height 100cm Height (folded) 114cm Height (folded) 124cm Width 42cm Width 41cm Width (folded) 42cm Width (folded) 41cm Depth 20cm Depth 22cm Weight 8.5Kg Weight 10.1Kg Max. Load 114Kg Max. Load 114Kg