



Retmex Rescue Sheets

For the rapid rescue of people from fires and hazardous locations, where time is critical. With a range of models for differing budgets and requirements.....



Retmex Compact



Retmex Comfort, with ergonomic carry handles



Retmex Bariatric



Retmex Nomex rescue sheet / fire blanket

The Retmex range of rescue sheets are designed for the quick and efficient rescue from a fire or from other hazardous situations.

Available in a range of models, from the budget-priced **Retmex Compact**, simple in design but effective in function, to the **Retmex Comfort** which includes additional features such as ergonomic carry handles, and the **Retmex Bariatric** which has been load-tested to evacuate a person up to 385 kg (60 stone). The flagship model, the fire-resistant **Retmex Nomex**, is designed specifically as a fire rescue sheet, also doubling up as a fire blanket.



Retmex Rescue Sheet Range



Features

- **Retmex Compact** – as its name infers, this is a compact basic rescue sheet which simply attaches to a rescuer’s waist using the integrated carabiner, available to immediately deploy as and when required. With hand holds at various points on the rescue sheet and a pocket at the foot-end to hold and support the person on the sheet, casualties can be quickly and easily extracted from dangerous or hazardous locations. The sheet wraps up in itself for easy storage / carrying.



- **Retmex Comfort** – based on the Retmex Compact, but including more handle positions and the inclusion of PVC-covered grip handles for improved comfort and grip. This version comes with a carry bag for storage.



- **Retmex Bariatric** – specially designed for moving larger-size casualties, it has a load capacity of 385 kg, dimensions of 2200 mm by 1200 mm, and total of 18 carry handles, allowing the safe movement of the casualty using a suitable number of rescuers.



Specifications

Retmex Compact p/n 1099876610	2050 mm	800 mm	250 kg / 39 st
Retmex Comfort p/n 10PVC76590	2050 mm	800 mm	250 kg / 39 st
Retmex Bariatric p/n 109960219	2200 mm	1020 mm	385 kg / 60.5 st
Retmex Nomex p/n 1099602001	2000 mm	780 mm	250 kg / 39 st

- **Retmex Nomex** – designed for fire rescue and other applications. Made from fire-resistant materials Nomex® and Kevlar, it doubles as both a rescue sheet and also a fire blanket for extinguishing flames. With its own storage/carry system consisting of an integrated strap & carabiner, the Retmex Nomex sheet can be attached to any fixing point for transporting and carrying and available to hand for quick deployment as required. Weighing only 1.5 kg it minimises the weight of equipment that needs to be carried.



Retmex Range - Features and Benefits

- Compact and convenient to carry, always to hand
- Quick to deploy, allowing immediate response and quick evacuation
- Lightweight but durable, offering a long service life
- Strong materials with reinforced strapping, offering a high load capacity
- Multiple carry handles ensuring good manual handling posture can be adopted
- Foot pocket for extra stability and support
- Easy to clean for enhanced hygiene standards
- Tested to DIN standards, conforming to EN 597-1, EN 597-2 (Nomex - EN531,ISO/FDIS11612,EN 469)
- Applicable to a range of sectors including Fire rescue, HART teams, Police, Rescue teams, estates, etc.
- Usable in a range of situations such as mass-casualty evac, building evac, terrorist attacks, extraction from riots or shootings

