

FERNO SAVER SAFE Evacuation Chair

Designed for use in an emergency, when you need to get somebody down and out of a building - this is the Evacuation Chair for the job.



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

Height Open	Height Folded	Depth	Width	Weight	Load Limit
1334 mm	1030 mm	914 mm	527 mm	11.2 kg	180 kg/28st/396lbs

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Saver Safe Evacuation Chair** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**

FERNO SAVER SAFE Evacuation Chair

This chair provides the total solution for any emergency evacuation situation.

Designed to meet the needs of 'The Regulatory Reform (Fire Safety) Order 2005' and 'The Equality Act 2010' this chair should ideally sit on each floor stairwell, on its bracket, with its dust cover - **ready for action when required, but hopefully never used.**

The **FERNO SAVER SAFE Evacuation Chair** is the total solution for any emergency evacuation situation. It glides down stairs easily and smoothly to ensure no heavy lifting or manual handling is required, although training of its use is essential for operators or responsible person(s).

This SAVER SAFE chair should be in multi-storey Public Buildings, Shopping Centres, Sport Arenas, Offices and Care Homes, in fact anywhere where peoples safety is the responsibility of others.

- **Safe easy transfer of patient** - seat cover held at all 3 sides of the chair and cushioned at the front ensures the transfer of person on and off the chair is safe and also comfortable for the occupant
- **Large weight capacity & physical size frame** - 28 stone safe working load and a large seat gives the chair a large occupant capacity capability
- **Improved operator control and comfort** – easy-grip rubber sleeves around the upper handle allow excellent grip when moving the chair and rubberised sleeve on the base cross-frame ensures a reliably strong and safe braking grip when preparing the chair for descent
- **Reassuring patient safety** – Adjustable head pad and patient retaining strap to enhance comfort and security
- **Better manoeuvrability** – wide rear kickstand adds strength and stability on the rear castors, preventing tipping of the chair and making horizontal movement of the occupant effortless
- **No marking** - light-coloured track ensures no marking of the stairs
- **Simple to clean and maintain** – removable seat cover with zips situated on the underside
- **Easy storage** – Lightweight, with compact profile it comes with wall bracket, dust cover and user guidance signage



Part Nos:

SAV-001

Saver Safe evacuation chair, including wall bracket, dust cover and user guide signage

FERNO - Saver Safe Evacuation Chair FPC005 10.20

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.