



**Wheel-Kit 71 W 32 for
Titan 32 Basket Stretcher**

WHEEL-KIT 71 W 32



Wheel-Kit 71 W 32 for Titan 32 Basket Stretcher

The wheel kit for the Titan 32 basket stretcher facilitates transport over longer distances.

With a few simple steps, the kit is quickly attached to the foot-end of the basket stretcher using 4 quick-release safety fasteners.

The 260 mm large wheels with all-terrain profile and pneumatic tyres ensure easy running even on uneven terrain.

Powder-coated aluminum frame.



GENERAL SPECIFICATIONS

	Length	Width	Height	Weight	Load Capacity
Wheel-Kit 71 W 32	700 mm	950 mm	480 mm	12 kg	400 kg