


FERNO EasyFix Vacuum Splints

Each splint moulds around the fracture for the ultimate in comfort and support



Leg Splint

Arm Splint

Full Leg Splint

FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Length	Width	Depth	Weight
1360 mm (Full Leg Splint)	730 mm (Full Leg Splint)	30 mm (Full Leg Splint)	1.55 kg
980 mm (Leg Splint)	610 mm (Leg Splint)	30 mm (Leg Splint)	1.15 kg
745 mm (Arm Splint)	390 mm (Arm Splint)	30 mm (Arm Splint)	0.75 kg
750 mm (set 1)	370 mm (set 1)	180 mm (set 1)	3.10 kg
750 mm (set 2)	370 mm (set 2)	180 mm (set 2)	4.65 kg
Pump Length	Pump Diameter	Hose Length	Pump Weight
241 mm	76 mm	610 mm	0.6 kg

FERNO – is a world leading innovator and manufacturer of patient transportation and rescue solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix Vacuum Splints** typifies our passion for top quality, highly functional design...

patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



FERNO EasyFix Vacuum Splints

Unbeatable stability for patients during transportation – with exceptional ease of application

EasyFix vacuum splints are as easy to apply as they are effective. The splint moulds around the individual fracture to deliver a high level of comfort while keeping the affected limb safely immobilised before treatment.

Ideal for use during transportation, EasyFix splints require a minimum of preparation and need not be removed for X-rays.

- **Improved design gives better linear rigidity** – using new channel granule chamber system
- **Exceptional stability** – granulates spread evenly to support the limb
- **Easy to use** – applied quickly with minimum preparation
- **Hygienic and fire retardant** – PVC-coated polyester fibre is easy to keep clean and disinfect
- **Versatile** – arm, leg and full leg splints available
- **Flexible and practical to use** – moulds easily around the injured limb, with no need to remove for X-rays
- **Multi-purpose** – full leg splint with extra VELCRO® strap can be used as a paediatric full body vacuum mattress



Full Leg Splint



Leg Splint



Arm Splint



Part Nos:

272101000	Full Leg Splint (136 cm)
272102000	Leg Splint (98 cm)
272103000	Arm Splint (74.5 cm)
272104000	Splint Set I – arm, leg, repair kit, pump and bag
272105000	Splint Set II – arm, leg, full leg, repair kit, pump and bag
272000701	Vacuum Pump

FERNO – EasyFix Vacuum Splints FPC006 10.15

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Ferno House, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.