

TRANSCEND™



The **Transcend™ Stair Chair** provides maximum safety when transporting a patient down stairs. Its improved options and features carry patient loads of up to 550 pounds. The chair and track system carry the weight, increasing patient safety and reducing operator fatigue and risk of injury.

- **5-position lift bar and handles** reduce strain of lifting and provide maximum control for operators of various heights
- **The PowerTraxx Conversion Kit** can retrofit powered tracks to the manual Transcend Stair Chair



EZ GLIDE® technology provide a comfortable ride on almost any surface with 60% less shock than competitive stair chairs.



Rear lift handles with improved ergonomics for lifting or carrying the chair.



IV clip is built into the lift bar to easily secure an IV bag to the chair.



Extending footrest provides added support and comfort for patients.

Transcend™ Stair Chair

Specifications



HEIGHT (MIN)	36 in 908 mm
HEIGHT (MAX)	60 in 1518 mm
WIDTH (OVERALL)	20 in 505 mm
WIDTH (SEAT)	16 in 408 mm
DEPTH (FRONT TO BACK)	30 in 760 mm
LENGTH (FOLDED)	36 in 912 mm
WIDTH (FOLDED)	20 in 507 mm
DEPTH (FOLDED)	7.9 in 200 mm
WEIGHT (MANUAL TRACK CHAIR)	32 lb 14.7 kg
WEIGHT (CARRY CHAIR)	28 lb 12.7 kg
LOAD LIMIT	550 lb 250 kg

Features & Benefits

- **All-metal frame construction** provides durability and long life
- **Smooth 5-inch front swivel wheels and extra tall 8-inch rear wheels** with integrated shock absorbing tires that navigate rough terrain and do not collect debris
- **One inch track-to-ground clearance** provides smooth rolling over carpet and rough surfaces, as well as easier descent down stairs
- **Extending footrest and lower track angle** increase patient sense of security
- **Patient seating surface** is contoured and has removable injection molded seat panels
- **Positive locking hinge mechanism** ensures chair is locked in place whether deployed for use or folded for carrying and storage



If you are seeking an **evacuation chair**, rather than a stair chair designed for use by EMS professionals, ask your FERNO sales representative about FERNO's line of evacuation chairs and solutions.

Accessories & Options

PATIENT HANDLING ACCESSORIES



Headrest
No. 0822023



Rear Lift Handles
No. 0822028



Restraint Set
No. 0822024



O2 Holder
Kit No. 0822053

TRANSCEND STAIR CHAIR OPTIONS



Wheel Locks



Extending Lift Bar



Footrest



Armrests



Tri-Spoke Wheels



Rear Lift Handles



Glide Bar