## **FERNO Paraguard**Carry Handles

For aiding the movement of casualties when using Paraguard stretcher







Designed to aid the transfer of casualties being rescued using a Paraguard stretcher, the Paraguard carry handles simply attach and lock on at each end of the stretcher, allowing two rescuers to safely move a casualty further distances with less effort, whilst also enhancing the manual handling methods deployed. Can also be used in conjuction with Paraguard shoulder harness, available as an optional add-on.

010772700 Paraguard Slot In Carry Handles

Weight 2.5 kg (2 off)
Load capacity 136 kgs / 21 stone / 300 lbs
Dimensions 54cm \* 36 cm \* 21 cm

FERNO - Paraguard carry handles FPC001 01.21



