

FERNO Transfer Board & Transfer Mattress

A simple, effective means of transferring or removing a body



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Product	Length	Width	Weight	Thickness
Transfer Board	1140 mm	460mm	2 kg	5mm
Transfer Mattress	1950mm	500mm	2.5 kg	35mm

FERNO – is a world leading innovator and manufacturer of patient transport and funeral equipment.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Transfer aids** typifies our passion for top quality, highly functional design... funeral & mortuary products that don't just do the job, they do it **BETTER**.


call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**
Mortuary

FERNO Transfer Board & **Transfer** Mattress

Ferno Transfer Board

A simple, effective means of transferring a body between adjoining surfaces including tables, beds or transport cots.



Transfer Mattress

Designed for dignity and safety of both the body and technician, the innovative features of this mattress include six carry handles to lift the patient onto any removal stretcher.

It has a substantial carry capacity of 220kg.



Part Nos:

010799200	Transfer Board
271008800	Transfer Mattress

FERNO - Funeral Transfer Aids FPC003 09.14

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**
Mortuary

Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.